

Lens availability charts



Free Form Progressives	Hoyalux iD MyStyle® 2	Hoyalux iD LifeStyle® 3*	Hoyalux Array® 2	Hoyalux Array® 2 Wrap	Amplitude® HD3/Mini
Minimum Fitting Height	Personalized to frame 14mm	18mm/14mm	21mm/19mm/17mm/14mm	21mm/19mm/17mm/14mm	18mm/14mm
Plastic 1.50	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)
Plastic 1.50 Polarized		-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.00) ●●●	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	
Plastic 1.50 Sensity® 2	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD
Plastic 1.50 Transitions®		-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X
Phoenix® 1.53	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)
Phoenix® 1.53 Coppertone®			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	
Phoenix® 1.53 Polarized		-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +3.00) ●●●	-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	
Phoenix® 1.53 Sensity® 2	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD
Phoenix® 1.53 Transitions®		-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X
BluTech 1.56 (Indoor/Outdoor)			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	
Poly 1.59	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +3.00)	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +3.00)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +3.00)
Clear Blue Filter 1.59			-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	
Poly 1.59 Polarized		-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +3.00) ●
Poly 1.59 Coppertone®			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●	
Poly 1.59 Sensity® 2	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +3.00) ●●/SD	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +3.00) ●●●/SD	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +3.00) ●●●/SD
Poly 1.59 Transitions®		-6.00 to +6.00 Cyl. to -4.00 (+1.00 to +3.00) ●●	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-6.00 to +4.00 Cyl. to -4.00 (+0.75 to +3.50) ●●/X/XPG	-6.00 to +6.00 Cyl. to -4.00 (+1.00 to +3.00) ●
Hi-Index 1.60	-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)
Hi-Index 1.60 Sensity® 2	-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD	-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD/SS		-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)/SD
Hi-Index 1.60 Transitions®			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-6.00 to +4.00 Cyl. to -4.00 (+0.75 to +3.50) ●●	
Hi-Index 1.67	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)
Hi-Index 1.67 Polarized		-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.00) ●●●	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	
Hi-Index 1.67 Sensity® 2	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD
Hi-Index 1.67 Transitions®		-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-6.00 to +4.00 Cyl. to -4.00 (+0.75 to +3.50) ●●/X/XPG	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50) ●●
Hi-Index 1.74	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)

X = Transitions® XTRActive® New Generation XPG = Transitions® XTRActive® Polarized™ Gray SF = Sensity® Fast Gray & Brown SD = Sensity® Dark SS = Sensity® Shine

*Not all iD LifeStyle 3 polarized lenses come with iD technology.

Lens availability charts



Free Form Progressives	Hoyalux Array®	Hoyalux Array® Wrap	Hoyalux Summit® ecp/cd BKS	Amplitude® BKS/Mini	HOYA GP Wide™
Minimum Fitting Height	21mm/19mm/17mm/14mm	21mm/19mm/17mm/14mm	19mm/14mm	19mm/14mm	18mm
Plastic 1.50	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)
Plastic 1.50 Polarized	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.50) ●●●
Plastic 1.50 Sensity® 2	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS
Plastic 1.50 Transitions®	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-5.00 to +4.00 Cyl. to -4.00 (+0.75 to +3.50) ●●/X/XPG	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG
Phoenix® 1.53	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)
Phoenix® 1.53 Coppertone®	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●
Phoenix® 1.53 Polarized	-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●
Phoenix® 1.53 Sensity® 2	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS
Phoenix® 1.53 Transitions®	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X	-5.00 to +4.00 Cyl. to -4.00 (+0.75 to +3.50) ●●/X	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X
BluTech 1.56 (Indoor/Outdoor)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)		
Poly 1.59	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)
Clear Blue Filter 1.59	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)
Poly 1.59 Polarized	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●●
Poly 1.59 Coppertone®	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●
Poly 1.59 Sensity® 2	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS
Poly 1.59 Transitions®	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-6.00 to +4.00 Cyl. to -4.00 (+0.75 to +3.50) ●●/X/XPG	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG
Hi-Index 1.60	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)
Hi-Index 1.60 Sensity® 2	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD/SS		-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD/SS
Hi-Index 1.60 Transitions®	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-6.00 to +4.00 Cyl. to -4.00 (+0.75 to +3.50) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●
Hi-Index 1.67	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50)
Hi-Index 1.67 Polarized	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●
Hi-Index 1.67 Sensity® 2	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS
Hi-Index 1.67 Transitions®	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-6.00 to +4.00 Cyl. to -4.00 (+0.75 to +3.50) ●●/X/XPG	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG
Hi-Index 1.74	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)

X = Transitions® XTRActive® New Generation XPG = Transitions® XTRActive® Polarized™ Gray SF = Sensity® Fast Gray & Brown SD = Sensity® Dark SS = Sensity® Shine

Lens availability charts



Premium SV, Bifocal & New Media Optics

iD Single Vision®

New MySV™
Single Vision
iQ™/SV

ST28 iQ™/ST28

iD Space™, iD
Screen™, iD Zoom™

Tact®
40/60 BKS

Sync III™

Minimum Fitting Height				Personalized to frame 18mm	18mm	
Plastic 1.50		-10.00 to +10.00 Cyl. to -6.00	-10.00 to +6.00 Cyl. to -4.00 (+0.75 to +3.50)	-8.00 to +5.00 Cyl. to -6.00 (+1.00 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00
Plastic 1.50 Polarized		-10.00 to +4.00 Cyl. to -6.00 ●●●				-10.00 to +6.00 Cyl. to -6.00 ●●●
Plastic 1.50 Sensity® 2		-10.00 to +4.00 Cyl. to -6.00 ●●●/SF/SD/SS		-8.00 to +5.00 Cyl. to -6.00 (+1.00 to +3.50) ●●●●	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 ●●●●/SF/SD/SS
Plastic 1.50 Transitions®		-8.00 to +9.00 Cyl. to -4.00 ●●/X/XP	-8.00 to +6.00 Cyl. to -4.00 (+1.00 to +3.00) ●●/X		-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XP	-10.00 to +6.00 Cyl. to -6.00 ●●/X/XP
Phoenix® 1.53		-13.00 to +9.00 Cyl. to -6.00	-8.00 to +6.00 Cyl. to -4.00 (+1.00 to +3.00) ●●/X	-8.00 to +5.00 Cyl. to -6.00 (+1.00 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00
Phoenix® 1.53 Coppertone®		-13.00 to +9.00 Cyl. to -6.00 ●●●				-10.00 to +6.00 Cyl. to -6.00 ●●
Phoenix® 1.53 Polarized		-12.00 to +6.00 Cyl. to -6.00 ●●●				-10.00 to +4.50 Cyl. to -6.00 ●●●
Phoenix® 1.53 Sensity® 2		-13.00 to +9.00 Cyl. to -6.00 ●●●●/SF/SD/SS		-8.00 to +5.00 Cyl. to -6.00 (+1.00 to +3.50) ●●●●/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 ●●●●/SF/SD/SS
Phoenix® 1.53 Transitions®		-13.00 to +9.00 Cyl. to -4.00 ●●/X	-9.00 to +7.00 Cyl. to -4.00 (+1.00 to +3.50) ●●		-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +6.00 Cyl. to -6.00 ●●/X
BluTech 1.56 (Indoor/Outdoor)		-8.00 to +7.00 Cyl. to -6.00			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00
Poly 1.59		-10.00 to +6.00 Cyl. to -6.00	-8.00 to +6.00 Cyl. to -4.00 (+1.00 to +3.00)		-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +6.00 Cyl. to -6.00
Clear Blue Filter 1.59		-10.00 to +6.00 Cyl. to -6.00			-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +6.00 Cyl. to -6.00
Poly 1.59 Polarized		-9.00 to +6.00 Cyl. to -6.00 ●●●				-13.00 to +6.00 Cyl. to -6.00 ●●●
Poly 1.59 Coppertone®		-9.00 to +6.00 Cyl. to -6.00 ●●●				-13.00 to +6.00 Cyl. to -6.00 ●●●
Poly 1.59 Sensity® 2		-9.00 to +6.00 Cyl. to -6.00 ●●●●/SF/SD/SS			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●●/SF/SD/SS	-13.00 to +6.00 Cyl. to -6.00 ●●●●/SF/SD/SS
Poly 1.59 Transitions®		-6.00 to +6.00 Cyl. to -4.00 ●●/X/XP			-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XP	-13.00 to +5.25 Cyl. to -6.00 ●●/X/XP
Hi-Index 1.60	-10.00 to +8.00 Cyl. to -6.00	-13.00 to +8.00 Cyl. to -6.00		-11.00 to +7.00 Cyl. to -6.00 (+1.00 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00
Hi-Index 1.60 Sensity® 2	-10.00 to +8.00 Cyl. to -6.00 ●●●●/SD/SS	-13.00 to +6.50 Cyl. to -6.00 ●●●●/SD/SS		-10.00 to +6.00 Cyl. to -6.00 (+1.00 to +3.50) ●●●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●●/SD/SS	-10.00 to +6.00 Cyl. to -6.00 ●●●●/SD/SS
Hi-Index 1.60 Transitions®		-10.00 to +8.00 Cyl. to -4.00 ●●			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +6.00 Cyl. to -6.00 ●●
Hi-Index 1.67	-15.00 to +10.00 Cyl. to -6.00	-15.00 to +8.00 Cyl. to -6.00		-13.00 to +7.00 Cyl. to -6.00 (+1.00 to +3.50)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +7.50 Cyl. to -6.00
Hi-Index 1.67 Polarized		-15.00 to +10.00 Cyl. to -6.00 ●●●				-13.00 to +6.00 Cyl. to -6.00 ●●●
Hi-Index 1.67 Sensity® 2	-15.00 to +10.00 Cyl. to -6.00 ●●●●/SD/SS	-15.00 to +8.00 Cyl. to -6.00 ●●●●/SF/SD/SS		-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●●/SF/SS	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●●●/SF/SD/SS	-13.00 to +7.50 Cyl. to -6.00 ●●●●/SF/SD/SS
Hi-Index 1.67 Transitions®		-10.00 to +8.00 Cyl. to -4.00 ●●/X/XP			-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XP	-13.00 to +7.50 Cyl. to -6.00 ●●/X/XP
Hi-Index 1.74	-15.00 to +10.00 Cyl. to -6.00 (1.70)	-15.00 to +10.00 Cyl. to -6.00 (1.70)				-13.00 to +8.00 Cyl. to -6.00 (1.74)

X = Transitions® XTRActive® New Generation XPG = Transitions® XTRActive® Polarized™ Gray SF = Sensity® Fast Gray & Brown SD = Sensity® Dark SS = Sensity® Shine

AR compatibility guide

Key	EX3+ & EX3				Recharge*				Super HiVision				HiVisionw/VP						
Warranty	2 Years/unlimited				2 Years/unlimited				2 Years/unlimited				2 Years/2 times						
Designs	Plastic	Plastic Sensity® 2	Plastic Transitions®	Plastic Polarized	Phoenix®	Phoenix® Sensity® 2	Phoenix® Transitions®	Phoenix® Polarized	Poly	Poly Sensity® 2	Poly Transitions®	Poly Polarized	Hi index	Hi index Sensity® 2	Hi index Transitions®	Hi index	Hi index Sensity® 2	Hi index Transitions®	Hi index
Refractive index	1.50	1.50	1.50	1.50	1.53	1.53	1.53	1.53	1.59	1.59	1.59	1.59	1.60	1.60	1.60	1.67	1.67	1.67	1.74
Hoyalux iD MyStyle® 2	■	■			■	■			■	■			■	■		■	■		■
Hoyalux iD LifeStyle® 3	■	■	■		■	■	■		■	■	■	■	■	■	■	■	■	■	■
Hoyalux Array® 2/W**	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
iD Space™																			
iD Screen™	■				■	■							■			■			
iD Zoom™																			
Tact® 40/60 BKS	■	■	■		■	■	■		■	■	■		■	■	■	■	■	■	
Sync III™	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
iD Single Vision™													■	■		■	■		■
MySV™ Single Vision IQ™	■	■	■	■	■	■	■	■	■	■	■	■		■		■	■	■	■
ST28 IQ™	■		■		■		■					■							
Refractive index	1.50	1.50	1.50	1.50	1.53	1.53	1.53	1.53	1.59	1.59	1.59	1.59	1.60	1.60	1.60	1.66	1.66	1.66	1.70
SV	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Aspheric-SV	■		■		■	■	■		■		■		■	■	■	■	■	■	■
ST28	■		■	■	■		■		■										
ST35	■			■					■										
RD22	■																		
RD24	■		■																
STT7X28	■		■	■					■										
STT8X35	■		■	■					■										

*Recharge, Super HiVision EX3+ and Super HiVision automatically include the View Protect Performance.
 **Array 2 Wrap 1.74 only available in PAV or PAR.