

iD WorkStyle® 3: The perfect lens for comfort and productivity

iD WorkStyle 3 shares the same technologies, 3D Binocular Vision™ and AdaptEase™ Technology, found in iD MyStyle® 3, designed to maximize comfort for specific work styles. iD WorkStyle 3 is the primary pair of glasses for patients who spend most of their day – whether at work or at play – in an office or up-close environment.

iD WorkStyle 3 is available in three designs:

- | | | |
|--|---|--|
| SPACE | SCREEN | ZOOM |
| <ul style="list-style-type: none"> • The primary focus is beyond the distance of a computer screen • Ideal for patients who are active in their office environment | <ul style="list-style-type: none"> • The primary focus is at an intermediate distance, like those working at a desk • Ideal for anyone working on dual monitors | <ul style="list-style-type: none"> • The primary focus is clarity at a short distance and widest possible near area • Ideal for hobbies and laptop work requiring high levels of concentration |



| Lens Material | Rx Range (CYL to -6.00) | Coating Compatibility and Treatment |
|-----------------|----------------------------------|--|
| 1.50 / 1.50 CBF | -8.00 to +5.00 (+0.75 to +3.50) | <p>AR: Compatible with all HOYA AR coatings. Super HiVision® Meiryō® EX4™, Super HiVision® EX3+™, Super HiVision® EX3™, Super HiVision®, HiVision ViewProtect®, HiVision®, Premium AR ViewProtect®, Premium AR, Hard Coat, Backside AR 1.74 (no Backside AR and no Hard Coat)</p> <p>Treatments: Sensity® 3, Sensity Dark, Sensity Shine, Sensity Colors</p> |
| 1.53 Phoenix® | -8.00 to +5.00 (+0.75 to +3.50) | |
| 1.59 / 1.59 CBF | -8.00 to +5.00 (+0.75 to +3.50) | |
| 1.60 | -11.00 to +7.00 (+0.75 to +3.50) | |
| 1.67 | -13.00 to +7.50 (+0.75 to +3.50) | |
| 1.74 | -13.00 to +8.00 (+0.75 to +3.50) | |

For more information contact your HOYA Business Consultant or visit www.hoyavision.us.



©2026 HOYA Corporation. All rights reserved. iD WorkStyle, iD MyStyle, Meiryō, Phoenix, Super HiVision, Super HiVision, Sensity, View Protect, Recharge, and HiVision are registered trademarks of HOYA Corporation. 3D Binocular Vision, AdaptEase, EX 3, EX3+, and EX4 are trademarks of HOYA Corporation. PN7620758 05.26