

Lens availability charts

Free Form Progressives	Hoyalux iD MyStyle® 2	Hoyalux iD LifeStyle® 3*	Hoyalux Array® 2	Hoyalux Array® 2 Wrap	Amplitude® HD3/Mini
Minimum Fitting Height	Personalized to frame 14mm	18mm/14mm	21mm/19mm/17mm/14mm	21mm/19mm/17mm/14mm	18mm/14mm
Plastic 1.50	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00)	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)
Clear Blue Filter 1.50			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	
Plastic 1.50 Polarized		-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	
Plastic 1.50 Sensity® 2	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS
Plastic 1.50 Transitions®			-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XP/PG	-5.00 to +4.00 Cyl. to -4.00 (+0.75 to +4.00) ●●/X/XP/PG	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●/X
Phoenix® 1.53	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50)	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00)	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)
Phoenix® 1.53 Coppertone®		-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +4.00) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	
Phoenix® 1.53 Polarized		-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +4.00) ●●	-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +4.00) ●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	
Phoenix® 1.53 Sensity® 2	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS
Phoenix® 1.53 Transitions®		-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●/X	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●/X	-5.00 to +4.00 Cyl. to -4.00 (+0.75 to +4.00) ●●/X	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●/X
BluTech 1.56 (Indoor/Outdoor)			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00)	
Poly 1.59	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +4.00)	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +4.00)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +4.00)
Clear Blue Filter 1.59			-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00)	
Poly 1.59 Polarized	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +4.00) ●
Poly 1.59 Coppertone®		-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●	
Poly 1.59 Sensity® 2	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +4.00) ●●●/SF/SD/SS	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +4.00) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-6.00 to +6.00 Cyl. to -6.00 (+1.00 to +4.00) ●●●/SF/SD/SS
Poly 1.59 Transitions®		-6.00 to +6.00 Cyl. to -4.00 (+1.00 to +4.00) ●●	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XP/PG	-6.00 to +4.00 Cyl. to -4.00 (+0.75 to +4.00) ●●/X/XP/PG	-6.00 to +6.00 Cyl. to -4.00 (+1.00 to +4.00) ●
Hi-Index 1.60	-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00)	-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)
Hi-Index 1.60 Polarized		-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00)	-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)
Hi-Index 1.60 Sensity® 2	-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD/SS		-11.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50) /SD
Hi-Index 1.60 Transitions®			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	-6.00 to +4.00 Cyl. to -4.00 (+0.75 to +4.00) ●●	
Hi-Index 1.67	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)
Clear Blue Filter 1.67			-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00)	
Hi-Index 1.67 Polarized		-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	
Hi-Index 1.67 Sensity® 2	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS
Hi-Index 1.67 Transitions®		-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XP/PG	-6.00 to +4.00 Cyl. to -4.00 (+0.75 to +4.00) ●●/X/XP/PG	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00) ●●/X
Hi-Index 1.74	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +3.50)

X = Transitions® XTRActive® New Generation XPG = Transitions® XTRActive® Polarized™ Gray SF = Sensity® Fast Gray & Brown SD = Sensity® Dark SS = Sensity® Shine
 *Not all iD LifeStyle 3 polarized lenses come with iD technology.

** Sensity Shine, Sensity 2 Blue/Green add powers are available from +0.75 to +3.50.

Lens availability charts

Free Form Progressives	Hoyalux Array®	Hoyalux Array® Wrap	Hoyalux Summit® ecp/cd BKS	Amplitude® BKS/Mini	HOYA GP Wide™
Minimum Fitting Height	21mm/19mm/17mm/14mm	21mm/19mm/17mm/14mm	18mm/14mm	19mm/14mm	18mm
Plastic 1.50	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)
Clear Blue Filter 1.50	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)
Plastic 1.50 Polarized	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●●
Plastic 1.50 Sensity® 2	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS
Plastic 1.50 Transitions®	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG	-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG
Phoenix® 1.53	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)
Phoenix® 1.53 Coppertone®	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●
Phoenix® 1.53 Polarized	-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +4.00) ●●	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +4.00) ●●	-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +4.00) ●●	-10.00 to +4.50 Cyl. to -6.00 (+0.75 to +4.00) ●●
Phoenix® 1.53 Sensity® 2	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS
Phoenix® 1.53 Transitions®	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●/X	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●/X	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●/X	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●/X
BluTech 1.56 (Indoor/Outdoor)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-5.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)		
Poly 1.59	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)
Clear Blue Filter 1.59	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)
Poly 1.59 Polarized	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●●
Poly 1.59 Coppertone®	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●
Poly 1.59 Sensity® 2	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS
Poly 1.59 Transitions®	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG	-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG
Hi-Index 1.60	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)
Hi-Index 1.60 Polarized	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)
Hi-Index 1.60 Sensity® 2	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD/SS		-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SD/SS
Hi-Index 1.60 Transitions®	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●
Hi-Index 1.67	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)
Clear Blue Filter 1.67	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)
Hi-Index 1.67 Polarized	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●	-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●
Hi-Index 1.67 Sensity® 2	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●/SF/SD/SS	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●●/SF/SD/SS
Hi-Index 1.67 Transitions®	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50) ●●/X/XPG	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XPG
Hi-Index 1.74	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)	-6.00 to +4.00 Cyl. to -6.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)

X = Transitions® XTRActive® New Generation XPG = Transitions® XTRActive® Polarized™ Gray SF = Sensity® Fast Gray & Brown SD = Sensity® Dark SS = Sensity® Shine
 ** Sensity Shine, Sensity 2 Blue/Green add powers are available from +0.75 to +3.50.

Lens availability charts

Premium SV, Bifocal & New Media Optics	iD Single Vision®	MySV™ Single Vision iQ™/SV	ST28 iQ™/ST28	iD Space™, iD Screen™, iD Zoom™	Tact® 40/60 BKS		Sync III™
					Personalized to frame 18mm	18mm	
Minimum Fitting Height							
Plastic 1.50		-10.00 to +10.00 Cyl. to -6.00	-10.00 to +6.00 Cyl. to -4.00 (+0.75 to +3.50)	-8.00 to +5.00 Cyl. to -4.00 (+1.00 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)		-10.00 to +6.00 Cyl. to -6.00
Clear Blue Filter 1.50		-10.00 to +10.00 Cyl. to -6.00	-10.00 to +6.00 Cyl. to -4.00 (+0.75 to +3.50)		-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)		-10.00 to +6.00 Cyl. to -6.00
Plastic 1.50 Polarized		-10.00 to +4.00 Cyl. to -6.00 ●●●					-10.00 to +6.00 Cyl. to -6.00 ●●●
Plastic 1.50 Sensity® 2		-10.00 to +4.00 Cyl. to -6.00 ●●●/SF/SD/SS		-8.00 to +5.00 Cyl. to -4.00 (+1.00 to +3.50) ●●●●	-8.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●●/SF/SD/SS		-8.00 to +6.00 Cyl. to -6.00 ●●●●/SF/SD/SS
Plastic 1.50 Transitions®		-8.00 to +9.00 Cyl. to -4.00 ●●/X/XP	-8.00 to +6.00 Cyl. to -4.00 (+1.00 to +3.50) ●●/X		-10.00 to +3.75 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XP		-10.00 to +6.00 Cyl. to -6.00 ●●/X/XP
Phoenix® 1.53		-13.00 to +9.00 Cyl. to -6.00	-8.00 to +6.00 Cyl. to -4.00 (+1.00 to +3.50) ●●/X	-8.00 to +5.00 Cyl. to -4.00 (+1.00 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)		-10.00 to +6.00 Cyl. to -6.00
Phoenix® 1.53 Coppertone®		-13.00 to +9.00 Cyl. to -6.00 ●●					-10.00 to +6.00 Cyl. to -6.00 ●●
Phoenix® 1.53 Polarized		-12.00 to +6.00 Cyl. to -6.00 ●●					-10.00 to +4.50 Cyl. to -6.00 ●●
Phoenix® 1.53 Sensity® 2		-13.00 to +9.00 Cyl. to -6.00 ●●●●/SF/SD/SS		-8.00 to +5.00 Cyl. to -4.00 (+1.00 to +3.50) ●●●●/SD/SS	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●●/SF/SD/SS		-10.00 to +6.00 Cyl. to -6.00 ●●●●/SF/SD/SS
Phoenix® 1.53 Transitions®		-13.00 to +9.00 Cyl. to -4.00 ●●/X	-9.00 to +7.00 Cyl. to -4.00 (+1.00 to +3.50) ●●		-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XP		-10.00 to +6.00 Cyl. to -6.00 ●●/X/XP
BluTech 1.56 (Indoor/Outdoor)		-8.00 to +7.00 Cyl. to -6.00			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)		-10.00 to +6.00 Cyl. to -6.00
Poly 1.59		-10.00 to +6.00 Cyl. to -6.00	-8.00 to +6.00 Cyl. to -4.00 (+1.00 to +3.50)		-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)		-13.00 to +6.00 Cyl. to -6.00
Clear Blue Filter 1.59		-10.00 to +6.00 Cyl. to -6.00	-8.00 to +6.00 Cyl. to -4.00 (+1.00 to +3.50)		-13.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)		-13.00 to +6.00 Cyl. to -6.00
Poly 1.59 Polarized		-9.00 to +6.00 Cyl. to -6.00 ●●●	-8.00 to +6.00 Cyl. to -4.00 (+1.00 to +3.50) ●●●				-13.00 to +6.00 Cyl. to -6.00 ●●●
Poly 1.59 Coppertone®		-9.00 to +6.00 Cyl. to -6.00 ●●●					-13.00 to +6.00 Cyl. to -6.00 ●●●
Poly 1.59 Sensity® 2		-9.00 to +6.00 Cyl. to -6.00 ●●●●/SF/SD/SS			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●●●/SF/SD/SS		-13.00 to +6.00 Cyl. to -6.00 ●●●●/SF/SD/SS
Poly 1.59 Transitions®		-6.00 to +6.00 Cyl. to -4.00 ●●/X/XP			-13.00 to +5.25 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XP		-13.00 to +5.25 Cyl. to -6.00 ●●/X/XP
Hi-Index 1.60	-10.00 to +8.00 Cyl. to -4.00	-13.00 to +8.00 Cyl. to -6.00	-10.00 to +8.00 Cyl. to -4.00 *Only ST28	-11.00 to +7.00 Cyl. to -4.00 (+1.00 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)		-10.00 to +6.00 Cyl. to -6.00
Hi-Index 1.60 Polarized	-10.00 to +8.00 Cyl. to -4.00	-13.00 to +8.00 Cyl. to -6.00	-10.00 to +8.00 Cyl. to -4.00 *Only ST28	-11.00 to +7.00 Cyl. to -4.00 (+1.00 to +3.50)	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00)		-10.00 to +6.00 Cyl. to -6.00
Hi-Index 1.60 Sensity® 2	-10.00 to +8.00 Cyl. to -4.00 ●●●●/SD/SS	-13.00 to +6.50 Cyl. to -6.00 ●●●●/SD/SS		-10.00 to +6.00 Cyl. to -4.00 (+1.00 to +3.50) ●●●●	-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +3.50) ●●●●/SD/SS		-10.00 to +6.00 Cyl. to -6.00 ●●●●/SD/SS
Hi-Index 1.60 Transitions®		-10.00 to +8.00 Cyl. to -4.00 ●●			-10.00 to +6.00 Cyl. to -6.00 (+0.75 to +4.00) ●●		-10.00 to +6.00 Cyl. to -6.00 ●●
Hi-Index 1.67	-15.00 to +10.00 Cyl. to -4.00	-15.00 to +8.00 Cyl. to -6.00		-13.00 to +7.00 Cyl. to -4.00 (+1.00 to +3.50)	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)		-13.00 to +7.50 Cyl. to -6.00
Clear Blue Filter 1.67		-15.00 to +8.00 Cyl. to -6.00			-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00)		-13.00 to +7.50 Cyl. to -6.00
Hi-Index 1.67 Polarized		-15.00 to +10.00 Cyl. to -6.00 ●●●	-15.00 to +10.00 Cyl. to -4.00 ●●●				-13.00 to +6.00 Cyl. to -6.00 ●●●
Hi-Index 1.67 Sensity® 2	-15.00 to +10.00 Cyl. to -4.00 ●●●●/SD/SS	-15.00 to +8.00 Cyl. to -6.00 ●●●●/SF/SD/SS		-13.00 to +7.50 Cyl. to -4.00 (+0.75 to +3.50) ●●●●/SF/SS	-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●●●/SF/SD/SS		-13.00 to +7.50 Cyl. to -6.00 ●●●●/SF/SD/SS
Hi-Index 1.67 Transitions®		-10.00 to +8.00 Cyl. to -4.00 ●●/X/XP			-13.00 to +7.50 Cyl. to -6.00 (+0.75 to +4.00) ●●/X/XP		-13.00 to +7.50 Cyl. to -6.00 ●●/X/XP
Hi-Index 1.74	-15.00 to +10.00 Cyl. to -4.00	-15.00 to +10.00 Cyl. to -6.00		-13.00 to +8.00 Cyl. to -4.00 (+0.75 to +3.50)	-13.00 to +8.00 Cyl. to -6.00 (+0.75 to +4.00)		-13.00 to +8.00 Cyl. to -6.00

X = Transitions® XTRActive® New Generation XPG = Transitions® XTRActive® Polarized™ Gray SF = Sensity® Fast Gray & Brown SD = Sensity® Dark SS = Sensity® Shine
 ** Sensity Shine, Sensity 2 Blue/Green add powers are available from +0.75 to +3.50.

AR compatibility guide

Key	New		Super HiVision® Meiryō™ EX4™		Super HiVision® EX3+™ & EX3™		Recharge™		Super HiVision®		HiVision® w/VP										
Warranty	2 Years/unlimited				2 Years/unlimited				2 Years/unlimited		2 Years/unlimited		2 Years/2 times								
Designs	Plastic / Plastic CBF	Plastic Sensity® 2	Plastic Transitions®	Plastic Polarized	Phoenix®	Phoenix® Sensity® 2	Phoenix® Transitions®	Phoenix® Polarized	Poly / Poly CBF	Poly Sensity® 2	Poly Transitions®	Poly Polarized	Hi-Index 1.60	Hi-Index 1.60 Sensity® 2	Hi-Index 1.60 Transitions®	Hi-Index 1.60 Polarized	Hi-Index 1.67 / Hi-Index 1.67 CBF	Hi-Index 1.67 Sensity® 2	Hi-Index 1.67 Transitions®	Hi-Index 1.67 Polarized	Hi-Index 1.74
Refractive index	1.50	1.50	1.50	1.50	1.53	1.53	1.53	1.53	1.59	1.59	1.59	1.59	1.60	1.60	1.60	1.60	1.67	1.67	1.67	1.67	1.74
Hoyalux iD MyStyle® 2	Plastic Only	■			■	■			Poly Only	■	■	■	■	■			1.67 Only	■			■
Hoyalux iD LifeStyle® 3	Plastic Only	■	■	■	■	■	■	■	Poly Only	■	■	■	■	■	■	■	1.67 Only	■	■	■	■
Hoyalux Array® 2/Wrap Hoyalux Array/Wrap	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Summit® ecp/ cd BKS	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Amplitude® BKS/Mini	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Amplitude® HD3/Mini	Plastic Only	■			■	■	■	■	Poly Only	■	■	■	■	■	■	■	1.67 Only	■	■	■	■
iD Space™ iD Screen™ iD Zoom™	Plastic Only	■			■	■			Poly Only	■			■	■			1.67 Only	■	■		■
Tact® 40/60 BKS	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Hoyalux GP Wide®	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Sync III™	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
iD Single Vision™													■	■			■	■			■
MySV™ Single Vision iQ™	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
SV Aspheric SV	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ST28 ST28 iQ™	■		■		■		■		■				■								■
ST35	■								■												
RD24	■																				
STT7X28	■								■												
STT8X35	■								■												

*Super HiVision Meiryō EX4, Super HiVision EX3+, Recharge, and Super HiVision automatically include the View Protect Performance.
 **Array 2 Wrap 1.74 only available in PAV or PAR.