

Hoyalux iD[®] progressive lenses

Experience



Visual



Harmony



Hoyalux iD MyStyle[®] 2 is the most personalized and technologically advanced design available today

HOYA understands that everyone sees differently. It is essential to not only take into account the prescription and required correction, but also the individual characteristics of the wearer. iD MyStyle 2 takes into account all the variables that are specific for each wearer.

Our patented Integrated Dual Surface™ (iD) design shapes lenses exclusively for the way your patients see. It's the first to shape a lens design on both the front and back surface. This results in more comfortable transitions between near and far distances, wider viewing zones, and a dramatic reduction in blurring in peripheral vision.

With Binocular Harmonization Technology™ (BHT), iD MyStyle 2 takes into account the differences between right and left prescriptions to balance the design performance. Thus, perfect and effortless focus, constant stability and excellent depth of vision are provided.



**Binocular
Harmonization
Technology**



Position of wear



iD design



Stable view



Infinite optics



**4 convenience
options**

4 convenience options

Customizing iD MyStyle 2 lens design just got easier – Choose from 4 convenient profiles to get you on your way to customizing the perfect lens for your patient.



Modern

Indoors/outdoors

- Anyone who prefers only one pair of lenses for convenience
- No clear preference for specific distance



Adventure

Outdoors

- Outdoor sports and activities
- Teacher/driver
- Walks dog daily
- Clear edge to edge distance



Detail

Indoors

- IT worker/mechanic/accountant
- Close work focus
- Puzzles, crafts, scrap booking, hobbyist



Stable

Indoors

- Experienced wearer
- Senior/sedentary/less active
- Eye gazer
- ST28 to PAL

Hoyalux iD LifeStyle® 3

Things couldn't get any clearer than with iD LifeStyle 3—HOYA's new and improved progressive design that comes with new variations matching your patient's lifestyle. iD LifeStyle 3 is the newest generation of the iD Integrated Dual Side™ (iD) progressive product family.

Design benefits

- Optimized visual acuity and a crisper image perception
- Noticeably more relaxed vision thanks to precise focusing
- More accurate correction

Available designs

- **iD LifeStyle 3 Indoor:** Designed for presbyopes who spend a lot of their time indoors, looking at digital devices or reading
- **iD LifeStyle 3 Urban:** Designed for presbyopes who are active, always on the go, and also enjoy reading books and using their digital devices
- **iD LifeStyle 3 Outdoor:** Designed for presbyopes who primarily wear their correction for distance use, and still require clear vision to see their phone and other near materials

Now available in polarized options in gray, brown and green

New



Traditional design



iD LifeStyle 3

Hoyalux Array® 2

ARRAY

2

HOYA's best-in-class technologies and quality continue to advance progressive lenses. Array 2 applies HOYA freeform backside surfacing technology to virtually any semi-finished single vision lens, in any material, allowing your practice to offer a greater selection of lens options. Array 2 expanded availability makes dispensing multiple pairs of glasses easy.

Available designs

- Hoyalux Array 2
VL, 11, 13, 15, 17
- Hoyalux Array 2 Wrap
VL, 11, 13, 15, 17

Benefits

- Vision stability
- Precision focus
- Ease of dispensing with variable length (VL)
- Effortless adaptation
- Higher customer satisfaction

www.hoyavision.us

©2022 HOYA Corporation. All rights reserved. Hoyalux iD MyStyle, Hoyalux iD LifeStyle, and Hoyalux Array are registered trademarks of HOYA Corporation. Integrated Dual Surface, and Binocular Harmonization Technology are trademarks of HOYA Corporation. 08.22

HOYA
FOR THE VISIONARIES

Simply
outstanding
progressive
lenses

