



| Range (D) | Effective Diameter (mm) | Nominal Base | TRUE Curve N=1.530 | Back Curve | Center Thickness (mm) | Edge Thickness (mm) |
|-----------|-------------------------|--------------|--------------------|------------|-----------------------|---------------------|
| NA | 75 | 6.00 | 5.83 | -4.57 | 8.5 | 6.7 |
| NA | 75 | 4.50 | 4.62 | -5.60 | 5.5 | 6.8 |
| NA | 75 | 4.00 | 4.11 | -6.00 | 4.5 | 7.1 |
| NA | 75 | 3.50 | 3.63 | -7.00 | 4.7 | 9.5 |
| NA | 75 | 2.00 | 2.13 | -8.09 | 4.5 | 13.7 |

Prescription Selection Chart



Fitting Height = 19mm min.