



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +6.00	70	6.00	5.99	-4.50	9.8	7.8
+2.25 to +4.00	75	5.00	5.07	-5.50	7.7	8.3
0.00 to +2.00	75	4.00	3.94	-5.50	5.4	7.6
-4.00 to -0.25	75	3.00	2.96	-5.50	6.2	9.8
-10.00 to -4.25	75	2.00	1.98	-7.00	7.6	14.5

Prescription Selection Chart



Fitting Height = 19mm min.