



| Range (D) | Effective Diameter (mm) | Nominal Base | TRUE Curve N=1.530 | Back Curve | Center Thickness (mm) | Edge Thickness (mm) |
|-----------------|-------------------------|--------------|--------------------|------------|-----------------------|---------------------|
| +4.25 to +6.00 | 70 | 6.50 | 6.52 | -4.50 | 9.0 | 6.0 |
| +2.25 to +4.00 | 75 | 5.50 | 5.38 | -5.50 | 7.8 | 7.4 |
| 0.00 to +2.00 | 75 | 4.50 | 4.44 | -5.50 | 6.2 | 7.2 |
| -4.00 to -0.25 | 75 | 3.50 | 3.47 | -7.00 | 7.1 | 11.8 |
| -10.00 to -4.25 | 75 | 2.00 | 1.99 | -7.00 | 8.6 | 14.7 |

Prescription Selection Chart



Fitting Height = 19mm min.