



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +6.00	70	6.50	6.52	-4.50	9.0	6.0
+2.25 to +4.00	75	5.50	5.38	-5.50	7.8	7.4
0.00 to +2.00	75	4.50	4.44	-5.50	6.2	7.2
-4.00 to -0.25	75	3.50	3.47	-7.00	7.1	11.8
-10.00 to -4.25	75	2.00	1.99	-7.00	8.6	14.7

Prescription Selection Chart



Fitting Height = 19mm min.