



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +6.00	75	6.00	5.99	-4.50	8.5	6.3
+2.25 to +4.00	75	5.00	5.07	-5.50	6.8	7.1
0.00 to +2.00	75	4.00	3.94	-5.50	5.0	7.1
-4.00 to -0.25	75	3.00	2.96	-5.50	5.3	8.8
-10.00 to -4.25	75	2.00	1.98	-5.50	7.2	11.2

Prescription Selection Chart



Fitting Height = 14mm min.