


Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +6.00	70	7.00	7.02	-5.50	10.1	8.0
+2.25 to +4.00	75	5.75	5.82	-5.50	9.1	8.0
0.00 to +2.00	75	4.50	4.58	-5.50	7.1	7.9
-4.00 to -0.25	75	3.50	3.55	-5.50	10.6	13.0
-10.00 to -4.25	75	2.00	2.12	-5.50	12.2	15.8

Prescription Selection Chart



Fitting Height = 14mm min.