



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +6.00	80	7.00	7.41	-5.50	14.8	10.0
+2.25 to +4.00	80	6.00	5.98	-5.50	12.0	9.8
-1.00 to +2.00	80	5.00	5.01	-5.50	10.4	9.7
-5.00 to -1.25	80	4.00	3.95	-5.50	10.2	11.5
-10.00 to -5.25	80	2.00	2.13	-5.50	10.8	15.1

Prescription Selection Chart



Fitting Height = 18mm min.