



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +5.00	75	6.50	6.68	-6.25	11.2	9.0
+2.25 to +4.00	75	5.50	5.72	-6.25	9.4	9.0
-1.00 to +2.00	75	4.50	4.44	-6.25	11.0	13.0
-5.00 to -1.25	75	3.50	3.72	-6.25	10.0	13.0
-6.00 to -5.25	75	2.00	2.12	-4.25	11.0	14.0

Prescription Selection Chart



Fitting Height = 18mm min.