



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +6.00	75	6.50	6.68	-5.50	14.1	17.2
+2.25 to +4.00	75	5.50	5.72	-5.50	13.7	15.8
-1.00 to +2.00	75	4.50	4.44	-5.50	11.9	13.3
-5.00 to -1.25	75	3.50	3.72	-5.50	17.2	16.6
-10.00 to -5.25	75	2.00	2.12	-5.50	15.4	14.6

Prescription Selection Chart



Fitting Height = 18mm min.