



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +6.00	75	6.00	5.99	-5.00	11.1	11.0
+2.25 to +4.00	75	5.00	5.00	-3.00	6.0	3.5
-1.00 to +2.00	75	4.25	4.35	-7.00	7.3	10.8
-4.00 to -1.25	75	3.50	3.69	-5.50	4.3	6.2
-10.00 to -4.25	75	2.00	2.16	-7.00	3.8	10.0

Prescription Selection Chart



Fitting Height = 18mm min.