



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+6.25 to +8.00	65	6.50	6.63	-3.50	9.9	6.8
+4.25 to +6.00	70	5.50	5.49	-5.50	9.8	9.9
+2.25 to +4.00	70	4.50	4.61	-5.50	8.2	9.3
-2.00 to +2.00	75	3.75	3.73	-7.00	7.3	11.8
-6.00 to -2.25	70	3.00	2.99	-5.50	3.6	6.4
-10.00 to -6.25	70	2.00	2.05	-7.00	4.5	10.1
-13.00 to -10.25	65	1.00	0.90	-9.50	6.8	15.9

Prescription Selection Chart



Fitting Height = 18mm min.