



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +5.00	70	6.50	6.75	-6.50	9.2	NA
+2.25 to +4.00	75	5.50	5.46	-6.50	9.1	NA
-1.00 to +2.00	75	4.25	4.27	-6.50	6.1	NA
-4.00 to -1.25	75	3.50	3.50	-8.00	5.5	NA
-6.00 to -4.25	75	2.00	2.02	-7.50	5.5	NA

Prescription Selection Chart



Right side  
Convex view  
75mm lens.



<< Design reference points >>

- ① Distance reference point.
- ② Near reference point.
- ③ Prism reference point.
- Thickness reference point.

Fitting Height = 18mm min.