


Glass (1.60)

HOYA


GP Wide Clear

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+2.25 to +4.00	75	5.00	5.08	-4.50	8.2	NA
-1.00 to +2.00	75	4.25	4.36	-5.25	6.0	NA
-4.00 to -1.25	75	3.50	3.53	-6.25	6.0	NA
-8.00 to -4.25	75	2.00	2.04	-6.50	6.0	NA

Prescription Selection Chart



Right side Convex view 75mm lens.



<< Design reference points >>

- ① Distance reference point.
- ② Near reference point.
- ③ Prism reference point.

Fitting Height = 18mm min.