



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +6.00	75	6.50	6.52	-5.00	9.9	7.7
+2.25 to +4.00	75	5.50	5.38	-5.50	7.6	7.6
-1.00 to +2.00	75	4.50	4.44	-5.50	6.1	7.4
-4.00 to -1.25	75	3.50	3.47	-6.75	5.7	9.7
-10.00 to -4.25	75	2.00	1.99	-6.75	8.9	15.0

Prescription Selection Chart



Fitting Height = 18mm min.