


Prescription Selection Chart

Power (D)	Diameter (mm)	Center Thickness* (mm)	Edge Thickness* (mm)
+5.00	65	4.9	1.0
+4.75	65	4.7	1.0
+4.50	65	4.5	1.0
+4.25	65	4.4	1.0
+4.00	65	4.1	1.0
+3.75	65	3.9	1.0
+3.50	65	3.7	1.0
+3.25	65	3.6	1.0
+3.00	65	3.4	1.0
+2.75	65	3.2	1.0
+2.50	65	3.0	1.0
+2.25	65	2.8	1.0
+2.00	70	2.8	1.0
+1.75	70	2.6	1.0
+1.50	70	2.4	1.0
+1.25	70	2.3	1.0
+1.00	70	2.1	1.2
+0.75	70	2.0	1.3
+0.50	70	1.9	1.4
+0.25	70	1.7	1.5
+0.00	75	1.6	1.6
-0.25	75	1.5	1.8
-0.50	75	1.5	2.0
-0.75	75	1.4	2.2
-1.00	75	1.3	2.4
-1.25	75	1.2	2.5
-1.50	75	1.2	2.8
-1.75	75	1.1	3.0
-2.00	75	1.0	3.2
-2.25	75	1.0	3.4
-2.50	75	1.0	3.7
-2.75	75	1.0	4.0
-3.00	75	1.0	4.3
-3.25	75	1.0	4.6
-3.50	75	1.0	4.8
-3.75	75	1.0	5.2
-4.00	75	1.0	5.4
-4.25	75	1.0	5.7
-4.50	75	1.0	5.9
-4.75	75	1.0	6.4
-5.00	75	1.0	6.6
-5.25	75	1.0	6.9
-5.50	75	1.0	7.2
-5.75	75	1.0	7.7
-6.00	75	1.0	7.9
-6.25	75	1.0	8.2
-6.50	75	1.0	8.4
-6.75	75	1.0	9.0
-7.00	75	1.0	9.3



*Ranges of up to +/-0.2mm for the same power but different cylinders. Check the Lens Description Standard (LDS) files for exact values.