



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.75 to +6.00	75	7.00	6.91	-5.50	10.3	8.4
+3.25 to +4.50	75	6.00	5.95	-5.50	8.6	7.9
+1.75 to +3.00	80	5.25	5.19	-5.50	9.9	9.5
-2.00 to +1.50	80	4.25	4.19	-5.50	8.2	9.3
-6.00 to -2.25	80	3.25	3.18	-5.50	10.9	13.6
-10.00 to -6.25	75	2.00	2.12	-5.50	14.6	18.6

Prescription Selection Chart



Fitting Height = 18 mm min.

