



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.75 to +6.00	75	6.75	6.80	-6.00	13.7	12.6
+3.25 to +4.50	75	5.75	5.85	-5.00	8.6	7.3
+1.75 to +3.00	75	5.00	4.91	-5.00	7.7	7.7
-2.00 to +1.50	75	4.00	3.96	-6.00	5.8	8.5
-6.00 to -2.25	75	3.00	2.97	-7.00	8.2	13.1
-10.00 to -6.25	75	2.00	1.99	-7.00	9.7	16.0

Prescription Selection Chart



Fitting Height = 18mm min.

