



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.75 to +6.00	77	6.25	6.32	-6.25	10.8	9.9
+3.25 to +4.50	77	5.50	5.51	-6.25	8.5	8.7
+1.75 to +3.00	77	4.75	4.66	-6.25	8.6	9.7
-2.00 to +1.50	77	4.00	4.02	-6.25	8.2	10.2
-6.00 to -2.25	77	3.00	2.91	-4.25	10.0	10.8

Prescription Selection Chart



Fitting Height = 18 mm min.

