

Prescription Selection Chart

Power (D)	Diameter (mm)	Center Thickness* (mm)	Edge Thickness* (mm)
+4.00	65	5.1	1.3
+3.75	65	4.8	1.3
+3.50	65	4.6	1.3
+3.25	65	4.3	1.3
+3.00	65	4.1	1.3
+2.75	65	3.8	1.3
+2.50	65	3.6	1.3
+2.25	65	3.4	1.3
+2.00	65	3.0	1.3
+1.75	65	2.7	1.3
+1.50	65	2.5	1.3
+1.25	65	2.2	1.3
+1.00	65	2.0	1.3
+0.75	65	2.0	1.5
+0.50	65	1.8	1.5
+0.25	65	1.5	1.5
+0.00	70	2.2	2.2
-0.25	70	2.0	2.3
-0.50	70	1.9	2.5
-0.75	70	1.8	2.6
-1.00	70	1.8	2.9
-1.25	70	1.6	3.0
-1.50	70	1.6	3.2
-1.75	70	1.4	3.4
-2.00	70	1.4	3.6
-2.25	70	1.4	3.9
-2.50	70	1.4	4.2
-2.75	70	1.4	4.4
-3.00	70	1.4	4.7
-3.25	70	1.3	5.0
-3.50	70	1.3	5.3
-3.75	70	1.3	5.6
-4.00	70	1.3	5.9
-4.25	70	1.3	6.1
-4.50	70	1.3	6.4
-4.75	70	1.3	6.8
-5.00	70	1.3	7.0
-5.25	70	1.3	7.3
-5.50	70	1.3	7.6
-5.75	70	1.3	7.8
-6.00	70	1.3	8.1

MINUS CYLINDER POWER

