



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.75 to +6.00	65	7.00	6.91	-5.50	12.6	11.1
+3.25 to +4.50	65	6.00	5.95	-5.50	7.0	6.4
+1.75 to +3.00	70	5.25	5.22	-5.50	6.8	7.0
-2.00 to +1.50	70	4.25	4.21	-5.50	6.8	8.2
-6.00 to -2.25	70	3.25	3.16	-5.50	9.5	11.4
-10.00 to -6.25	70	2.00	2.12	-5.50	14.7	18.0

Prescription Selection Chart



Fitting Height = 14mm min.

