

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+5.25 to +6.00	75	8.00	8.05	-5.76	10.8	7.1
+4.25 to +5.00	75	7.00	6.96	-6.25	8.2	7.1
+2.75 to +4.00	75	6.00	5.99	-5.76	8.1	7.8
+1.25 to +2.50	75	5.00	5.00	-6.25	6.2	8.1
-2.00 to +1.00	75	4.00	4.03	-5.76	7.1	9.6
-5.00 to -2.25	75	3.00	3.04	-5.76	6.0	9.9
-8.00 to -5.25	75	2.00	1.99	-5.76	8.2	13.5
-9.00 to -8.25	70	0.50	0.65	-5.76	8.5	14.7

Prescription Selection Chart

