

Plastic (1.50) **HOYA** SF Aspheric SV Clear Coated and Uncoated

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+7.75 to +8.00	65	8.00	8.82	-5.50	T1 9.0	6.4
+6.75 to +7.50	65	7.50	7.90	-5.50	T1 8.5	6.9
+5.75 to +6.50	70	6.50	6.92	-4.00	T1 9.2	6.4
+4.75 to +5.50	70	5.50	5.94	-4.00	T1 8.1	6.4
+3.75 to +4.50	70	5.00	5.20	-4.00	T1 7.3	6.4
+2.75 to +3.50	70	4.50	4.70	-4.00	T1 6.9	6.4
+1.75 to +2.50	70	4.00	4.20	-4.00	T1 6.4	6.4

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+0.75 to +1.50	70	3.50	3.68	-4.00	T1 6.0	6.4
-0.50 to +0.50	70	3.00	3.16	-4.00	T1 5.7	6.6
-1.50 to -0.75	70	2.50	2.64	-4.00	T1 5.3	6.7
-3.50 to -1.75	70	2.00	2.11	-5.50	T1 5.0	9.0
-5.50 to -3.75	70	0.75	0.79	-5.50	T1 6.0	11.3
-6.50 to -5.75	70	0.50	0.53	-4.00	T1 9.2	12.6
-8.50 to -6.75	70	0.40	0.53	-4.00	T1 13.1	16.6

Prescription Selection Chart

