

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+7.50 to +8.00	65	8.00	7.27	-1.75	8.0	2.9
+6.50 to +7.25	65	7.50	6.47	-1.75	7.4	3.2
+5.50 to +6.25	65	6.50	5.65	-1.75	6.6	3.2
+4.50 to +5.25	65	5.50	5.02	-1.75	5.8	3.0
+2.50 to +4.25	75	4.50	3.98	-4.50	4.9	6.1
+0.50 to +2.25	65	3.50	3.10	-5.50	6.1	8.8
-0.50 to +0.25	70	3.00	2.66	-2.50	3.5	3.5
-1.50 to -0.75	80	2.50	2.23	-7.00	6.1	13.9

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
-2.50 to -1.75	70	2.00	1.78	-5.50	5.0	9.4
-4.50 to -2.75	80	1.50	1.34	-4.50	4.0	8.5
-6.50 to -4.75	75	0.50	0.45	-4.50	6.0	11.0
-7.50 to -6.75	75	0.40	0.45	-7.00	4.0	12.8
-8.50 to -7.75	70	0.30	0.45	-7.00	8.9	16.6
-10.00 to -8.75	70	0.10	0.22	-7.00	12.3	20.3

Prescription Selection Chart

