

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+9.25 to +10.00	65	7.75	7.76	-3.50	8.4	4.7
+8.25 to +9.00	65	7.00	7.06	-3.50	7.7	4.8
+7.25 to +8.00	65	6.25	6.34	-3.50	7.0	4.8
+6.25 to +7.00	65	5.75	5.61	-3.50	6.3	4.8
+5.25 to +6.00	70	5.00	5.04	-3.25	6.5	5.0
+4.25 to +5.00	70	4.25	4.29	-3.25	5.7	5.1
+3.25 to +4.00	70	3.50	3.53	-3.25	4.8	5.0
+2.25 to +3.00	70	3.00	2.95	-3.25	4.4	5.1

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+1.25 to +2.00	75	2.50	2.56	-4.00	4.0	6.2
-1.00 to +1.00	75	2.00	1.98	-4.00	4.0	6.7
-3.00 to -1.25	75	1.25	1.38	-4.00	4.0	7.3
-5.00 to -3.25	75	0.80	0.79	-4.00	4.0	7.9
-8.00 to -5.25	75	0.50	0.59	-7.00	4.0	12.8
-10.00 to -8.25	70	0.30	0.59	-7.00	4.0	11.7
-15.00 to -10.25	65	0.10	0.59	-9.50	6.1	16.1

Prescription Selection Chart

