


Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+9.25 to +10.00	65	7.50	7.44	-3.50	T2 7.7	4.3
+9.25 to +10.00	65	7.50	7.44	-3.50	T3 11.3	7.9
+8.25 to +9.00	65	6.75	6.77	-3.50	T2 7.4	4.7
+8.25 to +9.00	65	6.75	6.77	-3.50	T3 10.6	7.9
+7.25 to +8.00	65	6.00	6.08	-3.50	T2 6.7	4.7
+7.25 to +8.00	65	6.00	6.08	-3.50	T3 9.8	7.8
+6.25 to +7.00	65	5.25	5.38	-3.50	T2 6.0	4.7
+6.25 to +7.00	65	5.25	5.38	-3.50	T3 10.1	8.8
+5.25 to +6.00	70	4.50	4.65	-3.25	T2 5.9	4.9
+5.25 to +6.00	70	4.50	4.65	-3.25	T3 9.5	8.5
+4.25 to +5.00	70	3.75	3.93	-3.25	T2 5.1	4.9
+4.25 to +5.00	70	3.75	3.93	-3.25	T3 8.6	8.4

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+3.25 to +4.00	70	3.00	3.19	-3.25	T2 4.3	4.9
+3.25 to +4.00	70	3.00	3.19	-4.50	T3 7.3	9.3
+2.25 to +3.00	70	2.50	2.64	-3.25	T2 4.0	5.1
+2.25 to +3.00	70	2.50	2.64	-3.25	T3 7.5	8.6
+1.25 to +2.00	75	2.00	2.08	-2.75	T2 5.3	6.3
+1.25 to +2.00	75	2.00	2.08	-4.75	T3 6.7	10.5
-0.25 to +1.00	75	1.75	1.70	-2.75	T2 4.1	5.5
-2.00 to -0.50	80	1.00	1.14	-1.50	T1 4.0	4.5
-2.00 to -0.50	80	1.00	1.14	-3.60	T3 6.6	10.3
-3.00 to -2.25	80	0.80	0.57	-1.50	T1 4.0	5.0
-4.00 to -3.25	80	0.70	0.57	-2.50	T1 4.0	6.7
-4.00 to -3.25	80	0.70	0.57	-6.00	T3 4.2	12.5

Prescription Selection Chart



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
-6.00 to -4.25	75	0.50	0.57	-3.25	T1 4.0	7.3
-6.00 to -4.25	75	0.50	0.57	-7.75	T3 4.0	14.2
-8.00 to -6.25	75	0.30	0.57	-5.00	T1 4.0	9.9
-8.00 to -6.25	75	0.30	0.57	-7.75	T3 6.4	16.6
-9.00 to -8.25	70	0.20	0.57	-7.00	T1 4.0	11.7
-10.00 to -9.25	70	0.10	0.57	-7.00	T1 4.0	11.8
-10.00 to -9.25	70	0.10	0.57	-8.75	T3 7.2	17.4
-15.00 to -10.25	65	0.00	0.57	-7.75	T1 4.0	11.6
-15.00 to -10.25	65	0.00	0.57	-9.50	T3 9.6	19.3