

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+7.25 to +8.00	65	8.00	7.87	-5.50	8.6	6.8
+6.25 to +7.00	65	7.00	6.98	-5.50	7.8	6.9
+5.25 to +6.00	70	6.00	6.05	-5.50	7.6	7.8
+4.25 to +5.00	70	5.00	5.12	-5.50	6.4	7.7
+3.25 to +4.00	70	4.75	4.66	-5.50	5.6	7.3
+2.25 to +3.00	70	4.25	4.18	-5.50	5.0	7.1
+1.25 to +2.00	75	3.75	3.70	-5.00	5.0	7.2
+0.25 to +1.00	75	3.25	3.21	-5.00	5.0	7.8

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
-1.00 to 0.00	75	3.00	2.97	-5.00	5.0	8.0
-2.00 to -1.25	75	2.25	2.23	-4.75	5.0	8.1
-3.00 to -2.25	75	1.75	1.74	-4.75	5.5	9.2
-4.00 to -3.25	75	1.25	1.24	-4.75	6.2	10.5
-5.00 to -4.25	75	1.00	0.99	-4.75	7.4	12.0
-6.00 to -5.25	75	0.70	0.75	-4.75	8.5	13.4
-7.00 to -6.25	70	0.60	0.75	-4.50	9.2	13.4
-8.00 to -7.25	70	0.50	0.75	-4.50	11.0	15.3
-10.00 to -8.25	70	0.40	0.75	-4.50	13.6	18.0

Prescription Selection Chart

