

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
NA	76	8.25	8.27	-6.00	12.3	9.0
NA	76	6.25	6.21	-6.00	8.7	9.0
NA	76	5.25	5.26	-6.00	8.7	10.5
NA	76	4.25	4.30	-4.00	10.6	10.5
NA	76	2.25	2.09	-4.00	10.6	13.6
NA	76	0.50	0.44	-4.00	10.6	16.0

Prescription Selection Chart

