

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +7.00	72	8.25	8.26	-6.00	9.4	6.0
+0.25 to +4.00	72	6.25	6.20	-6.00	8.5	8.2
-1.75 to 0.00	72	5.25	5.25	-6.00	8.5	9.5
-3.75 to -2.00	72	4.25	4.29	-4.00	8.5	8.1
-5.75 to -4.00	72	2.25	2.09	-4.00	9.0	11.4
-10.00 to -6.00	72	0.50	0.44	-4.00	9.0	13.5

Prescription Selection Chart

