

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+7.25 to +9.00	69	10.50	10.35	-6.33	14.8	9.5
+5.75 to +7.00	73	8.75	8.76	-6.33	12.0	8.4
+4.25 to +5.50	77	7.25	7.22	-6.33	13.6	12.0
+0.25 to +4.00	77	6.25	6.19	-6.33	13.6	13.7
-1.75 to 0.00	77	5.25	5.19	-4.22	11.9	10.3
-3.75 to -2.00	77	4.25	4.25	-4.22	14.8	14.8
-5.75 to -4.00	77	3.50	3.48	-2.03	18.8	16.7
-8.00 to -6.00	77	2.00	2.00	-2.03	16.7	16.8
-10.00 to -8.25	77	0.50	0.44	-4.00	10.6	15.8

Prescription Selection Chart

