

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+2.25 to +4.00	75	8.00	7.92	-5.50	15.5	11.7
-1.00 to +2.00	75	6.00	5.81	-5.50	13.1	12.8
-3.00 to -1.25	75	5.00	5.00	-5.50	12.0	12.9
-6.00 to -3.25	75	4.00	3.74	-5.50	13.7	16.4
-10.00 to -6.25	75	2.00	1.61	-2.50	18.5	19.6

Prescription Selection Chart

