

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +7.00	76	7.25	8.26	-6.00	12.80	11.20
+0.25 to +4.00	76	6.25	6.20	-6.00	12.80	12.90
-1.75 to 0.00	76	5.25	5.25	-6.00	11.10	9.60
-3.75 to -2.00	76	4.25	4.29	-4.00	14.0	14.00
-4.00 to -6.75	76	2.25	2.09	-4.00	15.90	16.00

Prescription Selection Chart

MINUS CYLINDER POWER

