

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+5.25 to +6.50	70	9.00	8.70	-7.00	17.0 T3	14.5
+3.75 to +5.00	70	8.00	7.81	-7.00	16.1 T3	14.9
+2.25 to +3.50	70	7.00	6.65	-7.00	14.1 T3	14.6
+0.75 to +2.00	75	6.00	5.85	-7.00	7.6 T2	9.4
-2.00 to +0.50	75	5.00	4.85	-7.00	4.5 T2	7.8
-4.00 to -2.25	75	4.00	3.88	-7.00	5.0 T2	9.7
-6.00 to -4.25	75	3.00	2.84	-6.00	4.6 T2	9.1
-8.00 to -6.25	70	2.00	1.93	-7.00	6.7 T2	13.0
-13.00 to -8.25	70	1.00	1.14	-7.00	8.7 T2	16.0

Prescription Selection Chart

