

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+8.75 to +8.00	70	8.50	8.41	-7.75	11.0	10.1
+4.25 to +6.50	70	7.75	7.70	-7.75	9.8	9.9
+2.25 to +4.00	70	6.25	6.24	-7.75	8.1	10.3
+0.25 to +2.00	70	4.75	4.74	-7.75	6.6	10.7
-1.50 to 0.00	75	4.00	3.97	-7.00	5.7	10.3
-4.00 to -1.75	75	3.25	3.37	-7.00	5.4	10.8
-8.00 to -4.25	75	2.50	2.40	-7.75	9.9	18.0
-15.00 to -8.25	70	0.50	0.30	-5.50	7.0	13.4

Prescription Selection Chart

