

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+6.75 to +8.00	70	8.50	8.41	-6.00	11.0	7.6
+4.25 to +6.50	75	7.75	7.70	-6.00	11.0	8.3
+2.25 to +4.00	75	6.25	6.24	-6.00	9.0	8.6
+0.25 to +2.00	75	4.75	4.74	-6.00	8.0	9.9
-1.50 to 0.00	75	4.00	3.97	-6.00	8.0	11.0
-4.00 to -1.75	75	3.25	3.37	-6.00	10.0	13.8
-8.00 to -4.25	75	2.50	2.40	-6.00	10.0	15.2
-15.00 to -8.25	70	0.50	0.30	-6.00	10.0	16.9

Prescription Selection Chart

