

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +6.50	74	8.25	8.26	-6.50	11.9	8.7
+2.75 to +4.00	74	6.25	6.25	-6.50	8.6	8.5
-3.00 to +2.00	74	4.25	4.24	-6.50	8.4	11.2
-8.00 to -3.25	74	2.25	2.26	-6.50	5.6	10.9

Prescription Selection Chart

