

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.25 to +6.00	76	8.00	8.05	-7.68	10.8	10.1
0.00 to +4.00	76	6.00	5.99	-6.25	8.1	8.5
-5.00 to -0.25	76	4.00	4.03	-5.76	7.1	9.8
-8.00 to -5.25	76	2.00	1.99	-5.76	8.2	13.8
-9.00 to -8.25	71	0.50	0.65	-5.76	8.5	14.5

Prescription Selection Chart

