

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+18.25 to +20.00	65 LC*	18.00	19.04	-7.75	14.5	12.8
+16.75 to +18.00	65 LC*	18.00	17.30	-7.75	16.7	13.5
+13.25 to +16.50	65 LC*	16.00	15.78	-7.75	16.5	14.1
+11.25 to +13.00	65 LC*	15.00	14.18	-7.75	17.0	14.6
+9.25 to +11.00	65	13.00	12.27	-7.75	20.0	13.6
+7.25 to +9.00	70	11.00	10.52	-7.00	12.5	6.9
+5.75 to +7.00	70	9.00	8.61	-5.50	13.7	9.5
+4.25 to +5.50	70	8.00	7.44	-8.75	7.9	9.8

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+2.25 to +4.00	70	7.00	6.53	-8.75	6.1	9.2
+0.25 to +2.00	70	6.00	5.46	-7.00	5.7	7.7
-1.75 to 0.00	70	5.00	4.47	-7.00	6.5	9.8
-3.75 to -2.00	70	4.00	3.49	-7.00	7.1	11.6
-5.75 to -4.00	70	3.00	2.53	-5.50	10.6	14.4
-10.00 to -6.00	70	2.00	1.52	-5.50	7.4	12.3
-13.00 to -10.25	65	1.00	0.54	-5.50	12.8	18.0

\*Lenticular Type

Prescription Selection Chart

