

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+5.25 to +6.00	77	8.00	8.05	-5.76	10.8	6.8
+4.25 to +5.00	77	7.00	6.98	-6.25	8.2	6.9
+2.75 to +4.00	77	6.00	5.99	-5.76	8.1	7.7
+1.25 to +2.50	77	5.00	5.00	-6.25	6.2	8.2
-2.00 to +1.00	77	4.00	4.03	-5.76	7.1	9.8
-5.00 to -2.25	77	3.00	3.08	-5.76	6.0	10.2
-8.00 to -5.25	77	2.00	1.99	-5.76	8.2	13.8

Prescription Selection Chart

