

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+5.25 to +6.00	76	8.00	8.05	-7.68	10.8	10.1
+4.25 to +5.00	76	7.00	6.97	-6.25	8.2	7.0
+2.75 to +4.00	76	6.00	5.99	-6.25	8.1	8.5
+1.25 to +2.50	76	5.00	5.00	-6.25	6.2	8.1
-2.00 to +1.00	76	4.00	4.06	-5.76	7.1	9.8
-5.00 to -2.25	76	3.00	3.04	-5.76	6.0	10.1
-8.00 to -5.25	76	2.00	2.02	-5.76	8.2	13.8
-9.00 to -8.25	71	0.50	0.71	-5.76	8.5	14.5

**Prescription Selection Chart**

