

| Range<br>(D)   | Effective<br>Diameter<br>(mm) | Nominal<br>Base | TRUE<br>Curve<br>N=1.530 | Back<br>Curve | Center<br>Thickness<br>(mm) | Edge<br>Thickness<br>(mm) |
|----------------|-------------------------------|-----------------|--------------------------|---------------|-----------------------------|---------------------------|
| +3.25 to +6.00 | 76                            | 8.25            | 8.26                     | -6.00         | 12.3                        | 9.0                       |
| 0.00 to +3.00  | 76                            | 6.25            | 6.21                     | -6.00         | 8.7                         | 9.0                       |
| -3.00 to -0.25 | 76                            | 4.25            | 4.30                     | -4.00         | 10.6                        | 10.5                      |
| -8.00 to -3.25 | 76                            | 2.25            | 2.09                     | -4.00         | 10.6                        | 13.5                      |

## Prescription Selection Chart

