

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+3.25 to +7.00	74	8.25	8.36	-6.25	11.5	8.5
-0.25 to +3.00	74	6.25	6.34	-6.25	8.0	8.0
-4.00 to -0.50	74	4.25	4.33	-6.25	12.0	14.5
-9.00 to -4.25	74	2.25	2.31	-4.25	12.5	14.0

Prescription Selection Chart

