

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+5.25 to +8.00	65	9.00	8.59	-5.50	13.2	9.7
+2.25 to +5.00	65	7.00	6.54	-5.50	11.0	10.0
-2.75 to +2.00	70	5.00	4.45	-7.00	9.4	12.7
-5.75 to -3.00	70	3.00	2.51	-7.00	12.4	18.0
-10.00 to -6.00	70	2.00	1.50	-7.00	11.2	18.0

Prescription Selection Chart

