



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.75 to +6.00	65	6.25	6.29	-5.50	10.5	9.6
+3.25 to +4.50	65	5.50	5.24	-4.50	6.4	5.1
+1.75 to +3.00	70	4.75	4.83	-5.50	5.5	6.1
-2.00 to +1.50	70	4.00	3.98	-5.50	5.5	7.2
-6.00 to -2.25	70	3.00	3.10	-7.00	5.5	9.5
-13.00 to -6.25	70	1.50	1.56	-7.75	5.5	12.3

Prescription Selection Chart



Fitting Height = 14mm min.