



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.75 to +6.00	65	6.00	5.91	-4.50	6.4	4.8
+3.25 to +4.50	65	5.25	5.22	-4.50	5.6	4.7
+1.75 to +3.00	65	4.50	4.50	-5.50	8.3	9.3
-2.00 to +1.50	65	3.75	3.77	-7.75	6.3	10.6
-6.00 to -2.25	65	3.00	3.03	-7.75	5.7	10.3
-10.00 to -6.25	65	2.00	2.09	-9.50	4.2	12.0
-13.00 to -10.25	65	1.00	0.95	-9.50	5.8	14.8

Prescription Selection Chart



Fitting Height = 14mm min.