



| Range (D) | Effective Diameter (mm) | Nominal Base | TRUE Curve N=1.530 | Back Curve | Center Thickness (mm) | Edge Thickness (mm) |
|-----------------|-------------------------|--------------|--------------------|------------|-----------------------|---------------------|
| +4.75 to +6.00 | 75 | 6.25 | 6.29 | -5.50 | 10.5 | 9.6 |
| +3.25 to +4.50 | 75 | 5.50 | 5.45 | -4.50 | 6.4 | 5.1 |
| +1.75 to +3.00 | 75 | 4.75 | 4.83 | -5.50 | 5.5 | 6.1 |
| -2.00 to +1.50 | 75 | 4.00 | 3.98 | -5.50 | 5.5 | 7.2 |
| -6.00 to -2.25 | 75 | 3.00 | 3.10 | -7.00 | 5.5 | 9.5 |
| -10.00 to -6.25 | 75 | 1.50 | 1.56 | -7.75 | 5.5 | 12.3 |

Prescription Selection Chart



Fitting Height = 14mm min.