



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.75 to +6.00	65	6.00	5.97	-5.50	9.1	10.3
+3.25 to +4.50	65	5.25	5.24	-5.50	9.0	8.9
+1.75 to +3.00	75	4.50	4.50	-7.00	8.9	12.3
-2.00 to +1.50	75	3.75	3.74	-5.50	5.1	7.4
-6.00 to -2.25	75	3.00	2.96	-7.00	6.7	11.6
-10.00 to -6.25	70	2.00	1.98	-7.00	9.0	15.8
-13.00 to -10.25	70	1.00	0.99	-7.00	11.0	17.7

Prescription Selection Chart



Fitting Height = 14mm min.